



Medical Product Catalogue

Our aim is to provide
innovative & humanized products and services



KNOW HONSUN ALL

Manufacturing Center

HONSUN's manufacturing center is located in Nantong. The city was given its name because it was the center of transportation in ancient China. And now, this city has become the foundation for the future of HONSUN.



The warehouse with an area of 4400 square meters

We built and run additional new warehouse with an area of 4,400 square meters. The warehouse is equipped with latest storage machinery and equipment.



OUR GROUP

HONSUN Group is a medical instrument enterprise group, which specializes in research, exploitation, production and distribution of health care product. Its headquarter locates in Shanghai, with the production bases of Wenzhou Yongqiang Seashore Industrial Zone and Nantong Economic Zone covering total area of 80,000 square meters.

We have been focusing on promoting the quality of people's life at healthcare by providing innovative and humanized products and services since founded in 1990. HONSUN is accredited with ISO13485 with products approved by CE & FDA according to AAMI-SP09: 2002, AAMI-SP10: 2002, EN1060-1/2/3. Reasonable price, excellent quality and perfect service are our strongpoint, accordingly the products are sold over 100 countries and regions in Europe, America, Oceania, Africa and Asia...

The powerful coexistence culture leads HONSUN to pursue opportune share with society, customers and staff, harmonious symbiosis, unremitting effort to improve human's health.



Heading Towards Industrial 4.0

Industrial automation is an important prerequisite for starting industry 4.0, Honsun, as the industry leader, is keeping pace with the trend of the times, We have been implementing automation on our core technology, and became the creator of advanced manufacture technology in the industry.



QUALITY CONTROL

In 1990, when we started our business, we decided that we would provide the best quality product to all of our clients. The story of HONSUN shows that we have done a great job in quality control, earning us a fantastic reputation.

We take quality control measures during every step of production. Because of the measures, we have obtained the CE certificate.



CONTENTS

ANEROID SPHYGMOMANOMETER ————— 04-11

HS-20 Aneroid Sphygmomanometer Series 05

HS-50 Aneroid Sphygmomanometer Series with Stethoscope 08

PALM TYPE SPHYGMOMANOMETER ————— 12-20

HS-201 Palm Type Sphygmomanometer Series 13

HS-GF-301 Sphygmomanometer For Emergency Use 20

WALL / DESK / STANDING ANEROID SPHYGMOMANOMETER ——— 21-28

HS-60 Wall/Desk Aneroid Sphygmomanometer Series 22

HS-70 Standing Aneroid Sphygmomanometer Series 24

STETHOSCOPE ————— 29-44

Cardiology Series 30

Clinical Series 37

Deluxe Series 39

Medical Series 41

NURSE AND STUDENT KIT ————— 45-46

Nurse Kit 45

Student Kit 46

ACCESSORIES ————— 47-52

HS-20 Aneroid Sphygmomanometer Series

As the advanced aneroid sphygmomanometer supplier in the world, we are glad to meet all of our customers' expectations with the most detailed features. Some features are not listed in this catalogue. However, we are experienced in developing new features and technology as requested.

HS-20F

Large Manometer Type

Technical Features	Model No.
ABS plastic case	1023
The dial diameter is Φ51.9mm	
Standard nylon cuff	132N-A
2-tube adult size pvc bladder	1222-A
Standard PVC bulb	144
Dual use air-release valve with spring	1528
Taper end valve	1554
Black drawstring bag	119N
Package Content	
Box size (mm): 108*65*168	
Carton size (mm) : 570*350*350	
Pcs/Ctn: 50	
Kgs/Ctn: 16.5/15	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	



HS-2000

Standard Type

Technical Features	Model No.
Special Aluminium alloy manometer	1068
Standard cotton cuff	132C-A
2-tube adult size PVC bladder	1222-A
Standard PVC bulb	144
Dual-use air-release valve with spring	1528
Taper end valve	1554
Vinyl zipper bag	113M
Package Content	
Box size (mm): 108*65*168	
Carton size (mm) : 570*350*350	
Pcs/Ctn: 50	
Kgs/Ctn: 20.5/19	
Additional adult cuff color options	
■ Black ■ Gray ■ Blue ■ Green	



HS-20A

Classic Type

Technical Features	Model No.
Special Aluminium alloy manometer	1068
Classic nylon cuff	133N-A
2-tube adult size PVC bladder	1222-A
Standard PVC bulb	144
Dual-use air-release valve with spring	1528
Taper end valve	1554
Vinyl zipper bag	113M
Package Content	
Box size (mm): 108*65*168	
Carton size (mm) : 570*350*350	
Pcs/Ctn: 50	
Kgs/Ctn: 22.5/21	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	



Plastic Storage Box
(optional)

HS-20D

DIY Type

Technical Features	Model No.
Special Aluminium alloy manometer	1068
Nylon cuff with D-ring	134N-A
2-tube adult size PVC bladder	1222-A
Standard PVC bulb	144
Dual-use air-release valve with spring	1528
Taper end valve	1554
Vinyl zipper bag	113M
Package Content	
Box size (mm): 108*65*168	
Carton size (mm) : 570*350*350	
Pcs/Ctn: 50	
Kgs/Ctn: 23.5/22	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	



Shockproof Ring PVC
(Optional)

HS-20C

Child Kit Type

Technical Features	Model No.
Special Aluminium alloy manometer	1068
Classic cartoon new born cotton cuff	133C-N
Classic cartoon child cotton cuff	133C-C
2-tube new born size PVC bladder	1222-N
2-tube child size PVC bladder	1222-C
Standard PVC bulb	144
Dual-use air-release valve with spring	1528
Taper end valve	1554
Cartoon zipper bag with ear	115C
Package Content	
Box size (mm): 110*75*185	
Carton size (mm) : 565*395*405	
Pcs/Ctn: 50	
Kgs/Ctn: 25/23.5	




HS-20Z

DIY Type

Technical Features	Model No.
Special Aluminium alloy manometer	1068-L
European Style nylon D-ring cuff	1341N-A
Standard PVC bulb	144
2-tube adult size PVC bladder	1222-A
Dual-use-air-release valve with spring	1528
Taper end valve	1554
Vlnyl zipper bag	113L
Package Content	
Box size (mm): 110*75*185	
Carton size (mm) : 565*395*405	
Pcs/Ctn: 50	
Kgs/Ctn: 17.5/16.5	
Additional adult cuff color options	
 Black 430U 2738U 296U 213U	





HS-50 Aneroid Sphygmomanometer Series with Stethoscope

As the advanced aneroid sphygmomanometer supplier in the world, we are glad to meet all of our customers' expectations with the most detailed features. Some features are not listed in this catalogue. However, we are experienced in developing new features and technology as requested.

HS-50A

Standard Sphygmomanometer Kit

Technical Features	Model No.
Special Aluminium alloy manometer	1068
Nylon cuff with D-ring	134N-A
2-tube adult size PVC bladder	1222-A
Standard PVC bulb	144
Dual-use air-release valve with spring	1528
Taper end valve	1554
Single head stethoscope	HS-30A
Large vinyl zipper bag	113L
Package Content	
Box size (mm): 110*75*185	
Carton size (mm) : 565*395*405	
Pcs/Ctn: 50	
Kgs/Ctn: 27.5/26	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	



HS-50B

Sphygmomanometer Kit With Dual Head Stethoscope

Technical Features	Model No.
Special Aluminium alloy manometer	1068
Classic nylon cuff	133N-A
2-tube adult size PVC bladder	1222-A
Standard PVC bulb	144
Dual-use air-release valve with spring	1528
Taper Metal net end valve (dust filter)	1555
Dual head stethoscope	HS-30B
Nylon carrying bag	116N
Package Content	
Box size (mm): 141*58*230	
Carton size (mm) : 480*440*325	
Pcs/Ctn: 30	
Kgs/Ctn: 18.5/17	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	



HS-50C

Home Use Sphygmomanometer Kit

Technical Features	Model No.
Special Aluminium alloy manometer	1068
Nylon cuff with D-ring (The stethoscope head is fixed to cuff)	1353N-A
2-tube adult size PVC bladder	1222-A
Standard PVC bulb	144
Dual-use air-release valve with spring	1528
Taper end valve	1554
Single head stethoscope	HS-30A
Large vinyl zipper bag	113L
Package Content	
Box size (mm): 141*58*230	
Carton size (mm) : 480*440*325	
Pcs/Ctn: 30	
Kgs/Ctn: 18.5/17	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	



HS-50D

Sphygmomanometer Kit With Rappaport Stethoscope

Technical Features	Model No.
Special Aluminium alloy manometer	1068
2-tube adult size PVC bladder	1222-A
Standard PVC bulb	144
Deflation valve: 'Push' for rapid deflation, 'screw' for precise deflation	1524
Taper Metal net end valve (dust filter)	1555
Colored rappaport stethoscope	HS-30C1
Name tag	173
Package Content	
Box size (mm): 141*58*230	
Carton size (mm) : 480*440*325	
Pcs/Ctn: 30	
Kgs/Ctn: 18.5/17	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	



HS-50E

Sphygmomanometer Kit With Large Manometer

Technical Features	Model No.
Large special Aluminium alloy manometer	1022
Nylon cuff with D-ring	134N-A
2-tube adult size PVC bladder	1222-A
Standard PVC bulb	144
Dual-use air-release valve with spring	1528
Taper end valve	1554
Single head stethoscope	HS-30A
Large vinyl zipper bag	113L
Package Content	
Box size (mm): 110*75*185	
Carton size (mm) : 565*395*405	
Pcs/Ctn: 50	
Kgs/Ctn: 29.5/28	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	




HS-50Z

Sphygmomanometer Kit With Deluxe single head Stethoscope

Technical Features	Model No.
Special Aluminium alloy manometer	1068-L
European Style nylon D-ring cuff	1341N-A
Standard PVC bulb	144
2-tube adult size PVC bladder	1222-A
Dual-use-air-release valve with spring	1528
Taper end valve	1554
Single head stethoscope	HS-30N
Nylon Carrying Case	116N
Package Content	
Box size (mm): 140*65*235	
Carton size (mm) : 450*350*500	
Pcs/Ctn: 30	
Kgs/Ctn: 18.5/16.5	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	





HS-201 Palm Type Sphygmomanometer Series

As the advanced aneroid sphygmomanometer supplier in the world, we are glad to meet all of our customers' expectations with the most detailed features. Some features are not listed in this catalogue. However, we are experienced in developing new features and technology as requested.

HS-201Z

Hollowed-out Design

Technical Features	Model No.
Unique Hollowed-out design case	1061
The thickness is 28.4mm	
The dial's diameter is 61.9mm	
Ultra-thin,one piece gauge, prefect design	
Newly designed air-release valve, beautiful appearance and good handle	1391
Nylon Bladder-less cuff	
Convenient disinfection	

Package Content
Box size (mm): 106*80*196.5
Carton size (mm) : 550*413*185
Pcs/Ctn: 20
Kgs/Ctn: 13/12kg




HS-201Y

Shock Resistant

Technical Features	Model No.
Unique 'shock resistant' case	1059
The dial's diameter is 55mm	
High-precision air- release valve	
Professional adult nylon cuff (anti-clamp)	135N-A
1-tube TPU bladder	1231-A
Flat bottom PVC bulb (special for 201Y)	
Black nylon bag	118N

Package Content
Box size (mm): 115*80*185
Carton size (mm) : 425*365*395
Pcs/Ctn: 30
Kgs/Ctn: 13.5/12

Additional adult cuff color options
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U



HS-201W

Self Adjustment

Technical Features	Model No.
ABS Plastic case and handle	1060
The dial's diameter is 62mm	
Standard adult Nylon cuff	132N-A
1-tube PVC bladder	1221-A
Standard PVC bulb (201S)	
Plastic end valve	1558
Special Vinyl Zipper case	113L
Package Content	
Box size (mm): 110*75*185	
Carton size (mm) : 565*395*405	
Pcs/Ctn: 50	
Kgs/Ctn: 16.5/15	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	



HS-201X

Good choice for women use

Technical Features	Model No.
Unique 'shock protected' case	1058
The dial's diameter is 55mm	
European style nylon cuff with D-ring	1341N-A
1-tube PVC bladder	1221-A
Special soft PVC bulb (special for 201X)	
Black nylon bag	118N
Package Content	
Box size (mm): 115*80*185	
Carton size (mm) : 425*365*395	
Pcs/Ctn: 30	
Kgs/Ctn: 13/11.5	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	



HS-201Q3

Technical Features	Model No.
Nickel plated metal case	1056
The dial's diameter is 59.5mm	
Two handling modes	
Zinc-alloy handle with Nickel plated	
Male & Female connector	
European Style nylon cuff with D-ring	1341N-A
1-tube PVC bladder	1221-A
Ribbed PVC bulb	147
Taper metal net end valve(dust filter)	1555
PVC Shockproof ring	1410
Black nylon bag	118N
Package Content	
Box size (mm): 115*80*185	
Carton size (mm) : 425*365*395	
Pcs/Ctn: 30	
Kgs/Ctn: 17.5/16	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	



HS-201Q1

Technical Features	Model No.
Nickel plated metal case	1051
The dial's diameter is 59.5mm	
Two handling modes	
Zinc-alloy handle with Nickel plated	
Deflation valve: 'Push' for rapid deflation, 'screw' for precise deflation	
Classic nylon cuff	133N-A
1-tube PVC bladder	1221-A
Ribbed PVC bulb	147
Taper metal net end valve(dust filter)	1555
Nylon carrying bag	116N
Package Content	
Box size (mm): 141*58*230	
Carton size (mm) : 480*440*325	
Pcs/Ctn: 30	
Kgs/Ctn: 20.5/19	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	



HS-201A1

Technical Features	Model No.
Chrome plated metal case	1041
The dial's diameter is 48mm	
Stainless steel handle	
Nylon cuff with D-ring	134N-A
1-tube PVC bladder	1221-A
Standard PVC bulb	144
Taper end valve	1554
Nylon zipper bag with ear	115N

Package Content

Box size (mm): 110*75*185
 Carton size (mm) : 565*395*405
 Pcs/Ctn: 50
 Kgs/Ctn: 25/23.5

Additional adult cuff color options

■ Black ■ 430U ■ 2738U ■ 296U ■ 213U



HS-201B2

Technical Features	Model No.
Chrome plated metal case	1057
The dial's diameter is 48 mm	
The professional configuration with twin tube technology	
Stainless steel handle with chrome plated	
Cotton Iron hook cuff for adult	1352C
2-tube adult latex bladder	1212-A
PVC bulb with end valve (special for 201B2)	
green nylon bag	118N

Package Content

Box size (mm): 115*80*185
 Carton size (mm) : 425*365*395
 Pcs/Ctn: 30
 Kgs/Ctn: 18/16.5

Options

Child hook cuff



HS-201J1

Technical Features	Model No.
ABS plastic upper case	1048
The dial's diameter is 59.5mm	
Stainless steel handle	
Cotton cuff with D-ring	134C-A
1-tube PVC bladder	1221-A
Standard PVC bulb	144
Taper end valve	1554
Cotton zipper bag with ear	115C

Package Content

Box size (mm): 110*75*185
 Carton size (mm) : 565*395*405
 Pcs/Ctn: 50
 Kgs/Ctn: 25/23.5

Additional adult cuff color options

■ Black ■ Gray ■ Blue ■ Green



HS-201K1

Technical Features	Model No.
ABS plastic upper case	1049
The dial's diameter is 60.2mm	
Two handling modes	
Large Zinc-alloy handle with black powder-coating	
Deflation valve: 'Push' for rapid deflation, 'screw' for precise deflation	
Nylon cuff with D-ring	134N-A
1-tube PVC bladder	1221-A
Ribbed PVC bulb	147
Taper metal net end valve(dust filter)	1555
Silicone shockproof ring	1412
Nylon carrying bag	116N

Package Content

Box size (mm): 141*58*230
 Carton size (mm) : 480*440*325
 Pcs/Ctn: 30
 Kgs/Ctn: 20.5/19

Additional adult cuff color options

■ Black ■ 430U ■ 2738U ■ 296U ■ 213U



HS-201T

Technical Features	Model No.
ABS plastic for case and handle	1053
The dial's diameter is 64.5mm	
Three handling modes	
European style Nylon cuff with D-ring	1341N-A
1-tube PVC bladder	1221-A
Flat bottom PVC bulb (special for 201T)	
Plastic end valve	1558
Nylon carrying bag	116N
Package Content	
Box size (mm): 141*58*230	
Carton size (mm) : 480*440*325	
Pcs/Ctn: 30	
Kgs/Ctn: 17/15.5	
Additional adult cuff color options	
■ Black ■ 430U ■ 2738U ■ 296U ■ 213U	



• Patent Protection

HS-201C1

Technical Features	Model No.
Chrome plated plastic case	1043
The dial's diameter is 60.2mm	
Stainless steel handle	
Cotton cuff with D-ring	134C-A
1-tube PVC bladder	1221-A
Standard PVC bulb	144
Taper end valve	1554
Large vinyl zipper bag	113L
Package Content	
Box size (mm): 110*75*185	
Carton size (mm) : 565*395*405	
Pcs/Ctn: 50	
Kgs/Ctn: 24/22.5	
Additional adult cuff color options	
■ Black ■ Gray ■ Blue ■ Green	



HS-201

Palm Type Sphygmomanometer Series (2-tube)

Package Content

Model	Box size (mm)	Carton size (mm)	Pcs/Ctn	Kgs/Ctn
HS-201A2	110*75*185	565*395*405	50	26.5/25
HS-201C2	110*75*185	565*395*405	50	25.5/24
HS-201Q2	141*58*230	480*440*325	30	22/20.5



HS-201A2



HS-201C2







HS-201Q2

HS-GF-301

Sphygmomanometer For Emergency Use

Technical Features	Model No.
Palm type sphygmomanometer	1043
Classic nylon for child/adult/large adult	133N-C/A/LA
1-tube PVC bladder for child	1221-C/A/LA
/adult/large adult	
Standard PVC bulb	144
Taper end valve	1554
Package Content	
Box size (mm): 265*235*55	
Carton size (mm) : 560*500*305	
Pcs/Ctn: 20	
Kgs/Ctn: 17.5/16	





HS-60/70 Wall/Desk/Standing Aneroid Sphygmomanometer Series


As the advanced aneroid sphygmomanometer supplier in the world, we are glad to meet all of our customers' expectations with the most detailed features. Some features are not listed in this catalogue. However, we are experienced in developing new features and technology as requested.



HS-60E

ABS Desk/ Wall Type

Technical Features	Model No.
Round shaped dial $\Phi 113$ mm	1074
Special precise mechanism	
The back plastic basket is a holder for cuff	
Nylon cuff with D-ring	134N-A
2-tube adult TPU bladder	1232-A
Standard latex bulb	142
Dual-use air-release valve with spring	1528
Standard end valve	1551
Spiral PVC tube (1.4m)	1413
Package Content	
Box size (mm): 138*138*175	
Carton size (mm) : 715*300*195	
Pcs/Ctn: 10	
Kgs/Ctn: 7.5/6	
Color Options	
 Black 430U 2738U 296U 213U	



HS-60A

ABS Desk/ Wall Type

Technical Features	Model No.
Square shaped dial 145mm×145mm	1070
Gradient adjustable for easy reading	
Standard cotton cuff	132C-A
2-tube adult latex bladder	1212-A
Standard latex bulb	142
Dual-use air-release valve with spring	1528
Standard end valve	1551
Spiral PVC tube (1.8m)	1413

Package Content

Box size (mm): 185*162*180
 Carton size (mm) : 820*390*205
 Pcs/Ctn: 10
 Kgs/Ctn: 10.5/9

Color Options

■ Black ■ Gray ■ Blue ■ Green



Fix on the wall

HS-60D

ABS Wall Type With ABS Basket

Technical Features	Model No.
Round shaped dial Φ147 mm	1073
Special precise mechanism	
180°direction adjustability for easy reading	
The back plastic basket is a holder of cuff	
Nylon cuff with D-ring	134N-A
2-tube adult latex bladder	1212-A
Standard latex bulb	142
Dual-use air-release valve with spring	1528
Standard end valve	1551
Spiral PVC tube (1.8m)	1413

Package Content

Box size (mm): 215*150*215
 Carton size (mm) : 795*445*240
 Pcs/Ctn: 10
 Kgs/Ctn: 12/10.5

Color Options

■ Black ■ 430U ■ 2738U ■ 296U ■ 213U



Fix on the wall

HS-70F

Stainless Steel Type With Metal Basket

Technical Features	Model No.
Square shaped dial 145mm×145mm	1080
Stainless steel castors & tube	
High stability	
360°swivel for easy reading	
Gradient adjustable for easy reading	
Adjustable height (90 cm to 133cm)	
Metal basket	
Cotton cuff with D-ring	134C-A
2-tube adult latex bladder	1212-A
Standard latex bulb	142
Dual-use air-release valve with spring	1528
Standard end valve	1551
Spiral PVC tube (1.8 m)	1413
Plastic connector	1567
Package Content	
Carton size (mm) : 750*220*205	
Pcs/Ctn: 1	
Kgs/Ctn: N/A	



HS-70E

Standing Type With Metal Basket

Technical Features	Model No.
Square shaped dial 145mm×145m	1079
Super-lower-gravity plastic castors	
High stability	
360°swivel for easy reading	
Gradient adjustable for easy reading	
Adjustable height (90 cm to 135 cm)	
Metal basket	
Cotton cuff with D-ring	134C-A
2-tube adult latex bladder	1212-A
Standard latex bulb	142
Dual-use air-release valve with spring	1528
Standard end valve	1551
Spiral PVC tube (1.8 m)	1413
Plastic connector	1567
Package Content	
Carton size (mm) : 800*320*167	
Pcs/Ctn: 1	
Kgs/Ctn: 4/3.2	



HS-70A

Standing Type With ABS Basket

Technical Features	Model No.
Square shaped dial 145mm×145mm	1075
Metal castors	
360°swivel for easy reading	
Gradient adjustable for easy reading	
5-wheel with brake	
Adjustable height (85 cm to 125 cm)	
Plastic basket	
Cotton cuff with D-ring	134C-A
2-tube adult latex bladder	1212-A
Standard latex bulb	142
Dual-use air-release valve with spring	1528
Standard end valve	1551
Spiral PVC tube (1.8 m)	1413
Plastic connector	1567

Package Content
Carton size (mm) : 575*495*160
Pcs/Ctn: 1
Kgs/Ctn: 5.5/4



HS-70D

Standing Type With ABS Basket

Technical Features	Model No.
Round shaped dial Φ147 mm	1078
Metal castors	
360°swiver for easy reading	
5-wheel with brake	
Adjustable height (95 cm to 135 cm)	
plastic basket	
Nylon cuff with D-ring	134N-A
2-tube adult latex bladder	1212-A
Standard latex bulb	142
Dual-use air-release valve with spring	1528
Standard end valve	1551
Spiral PVC tube (1.8m)	1413
Plastic connector	1567

Package Content
Carton size (mm) : 575*495*160
Pcs/Ctn: 1
Kgs/Ctn: 5.5/4



HS-70B

Standing Type With Metal Basket

Technical Features	Model No.
Square shaped dial 145mm×145mm	1076
Metal castors	
360°swivel for easy reading	
Gradient adjustable for easy reading	
5-wheel with brake	
Adjustable height (90 cm to 125 cm)	
Metal basket	
Cotton cuff with D-ring	134C-A
2-tube adult latex bladder	1212-A
Standard latex bulb	142
Dual-use air-release valve with spring	1528
Standard end valve	1551
Spiral PVC tube (1.8 m)	1413
Plastic connector	1567

Package Content
Carton size (mm) : 575*495*160
Pcs/Ctn: 1
Kgs/Ctn: 5.5/4





Stethoscope

Cardiology series: HS-106A , HS-107A and HS-108A are the featured stethoscopes in the Cardiology Series. Specially selected Materials and first-class craftsmanship ensure outstanding auscultation and comfortable use.

Clinical Series: Different spare parts can be used for multiple purpose.

Deluxe Series: Its reliable and clear acoustics offer an advantage to healthcare professionals when listening to the hearts, lungs and blood pressure sounds.

Medical Series: Economical design and good performance.



HS-106A

Master Cardiology Stethoscope

Technical Features


- Convenience is offered through HONSUN's "floating diaphragm" which allows you to hear low and high frequency sounds simply by applying light or heavy pressure to the chest-piece.
- Handcrafted solid stainless steel chestpiece. Excellent performance.
- The Non-chill rim and diaphragm is offer the comfort for your patient.
- One tube contains two parallel sound channels intended to enhance sound quality and minimize extraneous noise.

Specification

- Application: Cardiology/Excellent performance.
- Product weight/length: 215g/71cm.
- Chest piece material: stainless steel, Surface treatment: manual polish.
- Diaphragm: $\Phi 44\text{mm}$, Type: floating diaphragm.
- Binaural material: stainless steel, $\Phi 6\text{mm}$.
- Eartip material: silicone, Structure: soft sealing structure.
- Accessory box: 2 sets silicone ear tips.

Package Content


Box size (mm): 394*174*58
Carton size (mm) : 547*409*262
Pcs/Ctn: 12
Kgs/Ctn: 8/7



How to use Floating Diaphragm

To listen to low and high frequency sounds without removing and turning over the chestpiece, simply alternate between light and heavy pressure on the floating diaphragm.

Frequency response test



Additional tube color options

■ Black ■ 208C ■ 281C



HS-107A

Deluxe Cardiology Stethoscope

Technical Features


- The stethoscope adult side with HONSUN 'floating diaphragm'. As your requirement to adjustment frequency.
- Expect long-term durability with its corrosion, impact and scratch resistant stainless steel chest-piece.
- The Non-chill rim and diaphragm is offer the comfort for your patient.
- One tube contains two parallel sound channels intended to enhance sound quality and minimize extraneous noise.

Specification

- Application: Cardiology/High performance.
- Product weight/length: 203g/71cm.
- Chest piece material: stainless steel, Surface treatment: nature.
- Diaphragm: $\Phi 44\text{mm}/35\text{mm}$, Type: floating diaphragm/textured diaphragm.
- Binaural material: stainless steel, $\Phi 6\text{mm}$.
- Eartip material: Silicone, Structure: soft sealing structure.
- Accessory box: 2 sets silicone ear tips, Textured diaphragm 1 piece.

Package Content


Box size (mm): 394*174*58
Carton size (mm) : 547*409*262
Pcs/Ctn: 12
Kgs/Ctn: 8/7



High frequency
(child)



High frequency



Additional tube color options

■ Black ■ 208C ■ 281C



HS-108A

Deluxe Cardiology Stethoscope

Technical Features

- Classic Cardiology stethoscope is a unique, high-quality stethoscope specifically designed for use on Adult patients.
- A traditional bell on one side for lower frequency.
- The professional floating diaphragm on one side for higher frequency.
- Non-chill rim and diaphragm for comfort use.

Specification

- Application: Cardiology/High performance.
- Product weight/length: 220g/71cm.
- Chest piece material: stainless steel, Surface treatment: nature.
- Diaphragm: $\Phi 44\text{mm}$, Type: floating diaphragm.
- Binaural material: stainless steel, $\Phi 6\text{mm}$.
- Eartip material: Silicone, Structure soft sealing structure.
- Accessory box: 2 sets silicone ear tips.

Package Content


Box size (mm): 394*174*58
Carton size (mm) : 547*409*262
Pcs/Ctn: 12
Kgs/Ctn: 8.5/7.5



Low frequency



High frequency



Additional tube color options

■ Black ■ 208C ■ 281C

HS-101A

Special Zinc Alloy Cardiology Stethoscope

Specification

- Special Zinc alloy chestpiece plating. Sensitive $\Phi 48\text{mm}$ membrane for excellent acoustic performance.
- The concave design of chestpiece makes the sound collected twice.
- The two-in-one PVC tube system is great for sound conducting.
- Membrane with PVC anti-chill ring for avoiding cold shock.
- Silicone ear tip makes the auscultation clear and comfortable.
- Delicate design and compact size (Patent Protection).
- Comfortable and fixible binaural.
- Includes: Two pair of spare ear tips and a replaceable membrane.

Package Content

Box size (mm): 394*174*58

Carton size (mm) : 547*409*262


Pcs/Ctn: 12

Kgs/Ctn: 8.5/7.5

Additional tube color options

■ Black ■ 208C ■ 281C

Cardiology Series



HS-103A

Special Zinc Alloy Cardiology Stethoscope

Specification

- Silver dual head chestpiece for cardiology diagnosis, the $\Phi 48\text{mm}$ membrane is sensitive for sound collecting.
- The small, deep, cone-shaped chestpiece is used for auscultating of cacophony, as well as slight and low frequency sounds.
- The two-in-one PVC tube system is great for sound conducting.
- Considerately dual anti-chill ring design for avoiding cold shock.
- Silicone ear tip makes the auscultation clear and comfortable.
- Delicate design and compact size.
- Comfortable and fixible binaural.
- Includes: Two pair of spare ear tips and a replaceable membrane.

Package Content

Box size (mm): 394*174*58

Carton size (mm) : 547*409*262


Pcs/Ctn: 12

Kgs/Ctn: 10/9

Additional tube color options

■ Black ■ 208C ■ 281C

Cardiology Series





HS-30JZ

Shining Stethoscope

Specification

- The appearance design is attractive
- Shining tubing with new injection technology
- Chromed metal part is deluxe
- Stainless steel chestpiece with excellent performance
- comfortable with silicon eartip

Package Content


Box size (mm): 280*155*40
Carton size (mm) : 575*480*260
Pcs/Ctn: 36
Kgs/Ctn: 13.6/15.8

Additional tube color options

■ 493C ■ 452C ■ 7452C ■ 7717C



Pearlescent PVC Tube





HS-30NZ

Shining Stethoscope

Specification

- The appearance design is attractive
- Shining tubing with new injection technology
- Nickered metal part is textured
- Single-head design is more flexible
- comfortable with silicon eartip

Package Content

Box size (mm): 280*155*40
Carton size (mm) : 575*480*260
Pcs/Ctn: 36
Kgs/Ctn: 12.4/13.5

Additional tube color options

■ 493C ■ 452C ■ 7452C ■ 7717C



Pearlescent PVC Tube



HS-30J

Stainless Steel Cardiology Stethoscope

Cardiology Series



Specification

- Silver dual head chestpiece for adult, the $\Phi 48\text{mm}$ membrane is sensitive for sound collecting.
- The small bell-like chestpiece is used for auscultating slight and low frequency sounds.
- Considerately dual anti-chill ring design.
- Connecting spring is embedded, ear tips are specially designed to fit the structure of the ear.
- Stainless steel chestpiece, excellent sound conduction.
- Includes: A pair of Silicone ear trips

Package Content

Box size (mm): 252*132*38
Carton size (mm) : 690*530*230
Pcs/Ctn: 50
Kgs/Ctn: 15/13.5

Additional tube color options

■ Black ■ 1797C ■ 286C ■ 281C ■ 3288C ■ 298C
■ 212C ■ 1225C ■ 7C ■ 7677C



HS-30K

Stainless Steel Cardiology Stethoscope

Cardiology Series

Specification

- Silver dual head chestpiece for cardiology diagnosis, the $\Phi 48\text{mm}$ membrane is sensitive for sound collecting.
- The small, deep, cone-shaped chestpiece is used for auscultating of cacophony, as well as slight and low frequency sounds.
- Stainless steel chestpiece, excellent sound conduction.
- Considerately dual anti-chill ring design.
- 'A' shaped tube ensures the best sound conduction.
- Connecting spring is embedded, ear tips are specially designed for suiting the ear structure.
- Thick tube with excellent resistance to ambient interference.
- Includes: A pair of Silicone ear trips and a replaceable membrane.

Package Content

Box size (mm): 373*173*45
Carton size (mm) : 900*390*260
Pcs/Ctn: 25
Kgs/Ctn: 15/13.5

Additional tube color options

■ Black ■ 208C ■ 281C





HS-30C

Standard Sprague Rappaport Stethoscope

Specification

- Silver dual-membrane chestpiece: the big membrane with a diameter of 45.5mm is sensitive for sound collecting, the small membrane is applicable for child.
 - Zinc alloy chestpiece, excellent sound conduction, little acoustic attenuation.
 - Standard dual PVC tube, well-sealed dual sound conduction channel, good resistance to ambient interference, color is optional.
 - Attached transparent multi-pocket bag included.
 - Transparent PVC tube is optional
-

HS-30C1

Colored Sprague Rappaport Stethoscope

Specification

- Standard dual PVC tube, The color is optional. Attractive colored chestpiece with matching color tube.
 - Attached transparent multi-pocket bag included.
-

Package Content

Box size (mm): 300*130*40 Carton size (mm) : 670*620*230 Pcs/Ctn: 50 kgs/Ctn: 20/18.5


Additional tube color options

■ Black ■ 1797C ■ 021C ■ 3288C ■ 286C ■ 212C ■ 7677C ■ 7C

HS-30C2

Antimicrobial Sprague Rappaport Stethoscope

Clinical Series



Specification

- Antimicrobial function, patented design
- Antimicrobial silicone ear tips & tubing
- Efficiently kill Organism(s)
- Escherichia coli ATCC 8739
- Staphylococcus aureus ATCC 6538P
- More suitable for the hospital use
- More safety stethoscope
- TUV certificate
- Nylon zipper bag for stethoscope 120N

Package Content

Box size (mm): 286*164*50
Carton size (mm) : 840*306*275
Pcs/Ctn: 25
Kgs/Ctn: 15/13.5



HS-30C3

Clock Rappaport Stethoscope

Clinical Series

Specification

- Clock function, patented design.
- Silver single membrane chestpiece: the $\Phi 45.5\text{mm}$ membrane is sensitive for sound collecting.
- Zinc alloy chestpiece, good sound conduction, little acoustic attenuation.
- High quality rolled steel connecting spring, reasonable tightness.
- Accessory bag included.
- Transparent PVC bag packaging or gift box packaging.

Package Content

Box size (mm): 300*130*40
Carton size (mm) : 670*620*230
Pcs/Ctn: 50
Kgs/Ctn: 20/18.5

Additional tube color options

■ Black ■ 1797C ■ 021C ■ 3288C ■ 286C
■ 212C ■ 7677C ■ 7C



HS-30T1 *Clinical Series*

Standard Fetal Stethoscope

Specification

- Gynecological Silver chestpiece.
- Zinc alloy chestpiece, good sound conduction.
- Standard black dual PVC tube, well-sealed dual sound conduction channel, good resistance to ambient interference.
- High quality rolled steel connecting spring, reasonable tightness.

Package Content

Box size (mm): 252*132*38 Carton size (mm) : 690*530*230
Pcs/Ctn: 50 kgs/Ctn: 12.5/11



HS-30N *Deluxe Series*

Deluxe Single Head Stethoscope

Specification

- Silver single head chestpiece, ergonomic design for comfortable use. The diameter of chestpiece is 48mm, and it can catch both low frequency and high frequency sounds by pressing the back.
- Zinc alloy chestpiece with nickel-plated silver finish, excellent sound conduction, little acoustic attenuation.
- Connecting spring is embedded, ear tips are specially designed to fit the ear structure.
- Considerately anti-chill ring design.

Package Content

Box size (mm): 252*132*38 Carton size (mm) : 690*530*230
Pcs/Ctn: 50 kgs/Ctn: 14.5/14

Additional tube color options

Black 1797C 286C 281C 3288C 298C
212C 1225C 7C 7677C



HS-30L

Deluxe Dual Head Stethoscope

Specification

- Silver dual head chestpiece for adult, the $\Phi 47$ mm membrane is sensitive for sound collecting.
- The small bell-like chestpiece is used for auscultating slight and low frequency sounds.
- Connecting spring is embedded, ear tips are specially designed to fit the ear structure.
- Considerately dual anti-chill ring design.
- Aluminium alloy chestpiece, good sound conduction.

Package Content

Box size (mm): 252*132*38 Carton size (mm) : 690*530*230
Pcs/Ctn: 50 kgs/Ctn: 12/10.5

Additional tube color options

Black 1797C 286C 281C 3288C 298C
212C 1225C 7C 7677C





HS-30W

Sea Shell Stethoscope

Specification

- Delicate design and compact size (Patent Protection).
- Single head chestpiece for adult, the $\Phi 40\text{mm}$ membrane chestpiece is sensitive for sound collecting.
- Zinc alloy plus ABS with vacuum plating chestpiece, glossy and shiny.
- Considerately anti-chill ring design.
- Standard single PVC tube, color is optional.
- Fix with soft silicone ear tips.
- Includes: A pair of ear tips, a replaceable membrane and a anti-chill ring.

Package Content

Box size (mm): 260*130*38 Carton size (mm) : 675*275*415
Pcs/Ctn: 50 kgs/Ctn: 16/15

Additional tube color options

■ 648C ■ 2132C ■ 2095C ■ 186C ■ 1235C



HS-30Q

Cartoon Design Stethoscope

Specification

- Cartoon design stethoscope.
- The appearance is very kid-friendly. Reduce children's fear of doctor.
- Head special structure design, can change different head model.
- Considerately anti-chill ring design.
- Standard single PVC tube, color is optional.
- Fix with soft silicone ear tips.

Package Content

Box size (mm): 260*130*38
 Carton size (mm) : 675*275*415
 Pcs/Ctn: 50
 Kgs/Ctn: 14/13

Additional tube color options

Black 1797C 286C 281C 3288C 298C
 212C 1225C 7C 7677C



The head can change different model



Considerately anti-chill ring design



Soft silicone ear tips



Reduce children's fear of doctor



HS-30R1

Stethoscope For Neonate And Baby

Specification

- Silver dual head chestpiece for neonate and baby, The diameter of chestpiece is 23mm. Applicable for auscultating of the hearts and lungs of neonate and baby.
- Zinc alloy chestpiece, good sound conduction.
- Considerately dual anti-chill ring design.
- Standard single PVC tube, color is optional.

Package Content

Box size (mm): 195*100*28
 Carton size (mm) : 525*405*310
 Pcs/Ctn: 100
 Kgs/Ctn: 15.5/14

Additional tube color options

Black 1797C 286C 281C 3288C 298C
 212C 1225C 7C 7677C





HS-30A20

Sea Shell Stethoscope

Specification

- Single head chestpiece for adult, the $\phi 46\text{mm}$ membrane chestpiece is sensitive for sound collecting.
- Especially each chestpiece with different coloring coated design.
- Stainless steel chestpiece, good sound conduction.
- Considerately anti-chill ring design.
- Standard single PVC tube, color is optional.

Package Content

Box size (mm): 195*100*28
 Carton size (mm) : 525*405*310
 Pcs/Ctn: 100
 kgs/Ctn: 12/10.5

Additional tube color options

■ 648C ■ 2132C ■ 2095C ■ 186C ■ 1235C



🌾 Patent Protection: ZL202022316194.X 🌾



HS-30A / 30A2 / 30E

Single Head Stethoscope For Adult / Child



Specification

- Silver single-membrane chestpiece for adult, anodizing chestpiece, the $\Phi 47\text{mm}$ membrane is sensitive for sound collecting.
- Anodizing aluminum chestpiece for child, the $\Phi 38\text{mm}$ membrane is sensitive for sound collecting.
- Aluminium alloy chestpiece, good sound conduction.
- Standard single PVC tube, color is optional.

HS-30A1

Colored Single Head Stethoscope For Adult

Package Content

Model	Box size (mm)	Carton size (mm)	Pcs/Ctn	Kgs/Ctn
HS-30A	195*100*28	525*405*310	100	13.5/12
HS-30A1	195*100*28	525*405*310	100	13.5/12
HS-30A2	195*100*28	525*405*310	100	13/11.5
HS-30E	195*100*28	525*405*310	100	12.5/11

Additional tube color options

■ Black
 ■ 1797C
 ■ 286C
 ■ 281C
 ■ 3288C
 ■ 298C
 ■ 212C
 ■ 1225C
 ■ 7C
 ■ 7677C



HS-30B / 30B2 / 30F

Dual Head Stethoscope For Adult / Child



Specification

- Dual head chestpiece for adult. Anodizing aluminum chestpiece, the $\Phi 47\text{mm}$ membrane is sensitive for sound collecting. Applicable for auscultating of the heart and lungs.
- Dual head chestpiece for child. Anodizing aluminum chestpiece, the $\Phi 38\text{mm}$ membrane is sensitive for sound collecting. Applicable for auscultating of the hearts and lungs.
- The small bell-like chestpiece is used for auscultating slight and low frequency sounds.
- Aluminium alloy chestpiece, good sound conduction.
- Standard single PVC tube, color is optional.

HS-30B1 Colored Dual Head Stethoscope For Adult

Package Content

Model	Box size (mm)	Carton size (mm)	Pcs/Ctn	Kgs/Ctn
HS-30B	195*100*28	525*405*310	100	14.5/13
HS-30B1	195*100*28	525*405*310	100	14.5/13
HS-30B2	195*100*28	525*405*310	100	14.5/13
HS-30F	195*100*28	525*405*310	100	13.5/12

Additional tube color options

■ Black ■ 1797C ■ 286C ■ 281C ■ 3288C ■ 298C ■ 212C ■ 1225C ■ 7C ■ 7677C

SCK-001

Nurse Kit

Medical kit solution

Specification

Sturdy and tough First Aid Kit Bag with Multiple Compartments, Portable and Lightweight.
We are a professional manufacturer of medical tools, providing one-stop set matching services. Including medical student sets, ENT doctor sets, nurse sets, first aid sets, etc. You can configure different medical diagnostic tools according to your needs, such as pen lights, bandage scissors, Taylor hammers, multifunctional percussion hammers, tuning forks, etc. Or professional medical equipment, such as stethoscopes, oximeters, forehead thermometers, thermometers and other medical equipment, can be provided at one time.

Package Content

Box: 29*13.5*6.5/19.44kg Carton: 57*44*31/20.5kg Pcs/Ctn: 24

Color options



SCK-002

Student Kit

Medical kit solution

Specification

Sturdy and tough First Aid Kit Bag with Multiple Compartments, Portable and Lightweight.
We are a professional manufacturer of medical tools, providing one-stop set matching services. Including medical student sets, ENT doctor sets, nurse sets, first aid sets, etc. You can configure different medical diagnostic tools according to your needs, such as pen lights, bandage scissors, Taylor hammers, multifunctional percussion hammers, tuning forks, etc. Or professional medical equipment, such as stethoscopes, oximeters, forehead thermometers, thermometers and other medical equipment, can be provided at one time.

Package Content

Box: 29*13.5*6.5mm /14.08kg Carton: 61*41.5*41.5/16.5kg Pcs/Ctn: 32

Color options



Accessories



Specification	Model No.	Surface Treatment	Model No.
Special Aluminum alloy case	1068	Powder-coating finished	G
1068 Packing Bill	Box size (mm)	Carton size (mm)	Pcs/Ctn Kgs/Ctn
Foam box dimension	108*75*55	455*400*300	100 13.5/12.5



Specification	Size (mm)	Model No.
Special Vinyl Zipper Case (M)	160*100*55	113M
Special Vinyl Zipper Case (L)	170*100*65	113L



Specification	Size (mm)	Model No.
Cotton Zipper Case	170*100*65	115C
Nylon Zipper Case	170*100*65	115N



Specification	Size (mm)	Model No.
Cotton Carrying Case	235*125*50	116C
Nylon Carrying Case	235*125*50	116N
PVC-coated Carrying Case	235*125*50	116P



Specification	Size (mm)	Model No.
Nylon bag	190*110*70	118N



Specification	Size (mm)	Model No.
Black drawstring bag	150*200	119N



Specification	Size (mm)	Model No.
Nylon zipper bag for stethoscope	270*140	120N



Specification	Size (mm)	Model No.
Orange carrying bag	445*235*50	LD-FB-005



Specification	Size (mm)	Model No.
Plastic Storage Box	175*105*71	179

Accessories

Bladders



1211 1-tube




1212 2-tube

Product Name	Specification	Size (mm)	Model No.
Latex Bladder	New born A	40*80	1211/1212-NA
	New born B	30*90	1211/1212-NB
	Infant A	60*110	1211/1212-IA
	Infant B	50*190	1211/1212-IB
	Child	80*150	1211/1212-C
	Adult	120*220	1211/1212-A
	Large adult (Obese)	145*315	1211/1212-LA
	Thigh	180*370	1211/1212-T
	1212-A: 540*410*230 (mm) ; 200 PCS/CTN ; 24/23kg		
	1211-A: 540*410*230 (mm) ; 200 PCS/CTN ; 19.5/18.5kg		



1221 1-tube



1222 2-tube

Product Name	Specification	Size (mm)	Model No.
PVC Bladder	New born A	40*80	1221/1222-NA
	New born B	30*90	1221/1222-NB
	Infant A	60*120	1221/1222-IA
	Infant B	40*190	1221/1222-IB
	Child	80*150	1221/1222-C
	Adult	120*220	1221/1222-A
	Large adult (Obese)	145*315	1221/1222-LA
	Thigh	180*370	1221/1222-T
	1222-A: 540*410*230 (mm) ; 200 PCS/CTN ; 18/17kg		
	1221-A: 540*410*230 (mm) ; 200 PCS/CTN ; 13/12kg		



1231 1-tube




1232 2-tube

Product Name	Specification	Size (mm)	Model No.
TPU Bladder	New born A	40*80	1231/1232-NA
	New born B	30*90	1231/1232-NB
	Infant A	60*120	1231/1232-IA
	Infant B	40*190	1231/1232-IB
	Child	80*150	1231/1232-C
	Adult	120*220	1231/1232-A
	Large adult (Obese)	145*315	1231/1232-LA
	Thigh	180*370	1231/1232-T



1241 1-tube



1242 2-tube

Product Name	Specification	Size (mm)	Model No.
PVC Bladder	Infant	110*60	1241/1242-I
	Child	150*85	1241/1242-C
	Adult	225*120	1241/1242-A
	Large adult (Obese)	315*145	1241/1242-LA
	Thigh	370*185	1241/1242-T



Accessories

Cuff

Cuff		Material - C: Cotton		N: Nylon	
	Model No.	Specification	Size (mm)	Arm Range (mm)	
	132C/N-N	New born	250*45	---	
	132C/N-I	Infant	380*70	130 ~ 190	
	132C/N-C	Child	450*100	184 ~ 267	
	132C/N-A	Adult	500*140	254 ~ 406	
	132C/N-LA	Large adult (Obese)	600*170	333 ~ 510	
	132C/N-T	Thigh	780*215	406 ~ 660	
132C-A: 540*410*230 (mm) ; 200 PCS/CTN ; 11/10kg 132N-A: 540*410*230 (mm) ; 200 PCS/CTN ; 7/6kg					
	Model No.	Specification	Size (mm)	Arm Range (mm)	
	133C/N-N	New born	185*55	---	
	133C/N-I	Infant	255*75	130 ~ 190	
	133C/N-C	Child	345*110	184 ~ 267	
	133C/N-A	Adult	540*145	254 ~ 406	
	133C/N-LA	Large adult (Obese)	620*175	333 ~ 510	
	133C/N-T	Thigh	780*220	406 ~ 660	
133C-A: 540*410*230 (mm) ; 200 PCS/CTN ; 10/9kg 133N-A: 540*410*230 (mm) ; 200 PCS/CTN ; 8/7kg					
Special Cuff	Product Name	Specification	Size (mm)	Model No.	
	D-ring Cuff	Adult	500*140	134C/N-A	
		Large Adult (Obese)	600*170	134C/N-A	
134C-A: 540*410*230 (mm) ; 200 PCS/CTN; 14.5/13.5kg 134N-A: 540*410*230 (mm) ; 200PCS/CTN; 13.5/12.5kg					
	D-ring Cuff	European Style	500*150	1341C/N-A	
		1341N-A: 540*410*230 (mm) ; 200 PCS/CTN; 14/13kg			
	Nylon Bladderless Cuff	New born	250*50	1391/1392-N	
		Infant	215*55	1391/1392-I	
		Child	270*75	1391/1392-C	
		Adult	520*140	1391/1392-A	
		Large adult (Obese)	610*170	1391/1392-LA	
	Professional Adult Nylon Cuff	Adult	500*135	135N/A	

Accessories

Bulb			
Model No.	144	144A	142
			
	Standard PVC Bulb	Standard PVC Bulb With End Valve	Standard Latex Bulb

End Valve					
Model No.	1551	1552	1554	1555	1558
					
	Standard End Valve	Metal Net End Valve (Dust Filter)	Taper End Valve	Taper Metal Net End Valve (Dust Filter)	Insert Plastic End Valve

Air-release Valve			
Model No.	1523	1524	1528
			
Deluxe Air-release Valve	Push-button inching Valve (Precise Mini-deflation)	Dual-use Air-release Valve with Spring	

Connector			
Model No.	1565	1566	1567
			
	Male Connector	Female Connector	Plastic Connector

Accessories

Others

Model No.	1415	1410	1411	1412	1413
					
	Shockproof Ring PVC	PVC Shockproof Ring	Silicone Shockproof Ring (Small)	Silicone Shockproof Ring (Large)	PVC Coil Tube (1.4m/1.8M/2.0M/2.5M)

Accessory

Model No.	1713	1714	1715	1716
				
	"A" -shaped Tube with Stainless Steel Binaural	"Y" -shaped Tube with Stainless Steel Binaural	"Y" -shaped Tube with Copper Binaural	Two-in-one Tube System with Anatomically Shaped Binaural For HS-101A

Accessory

Model No.	144	144A	142	142
				
	Two-in-one Tube System with Anatomically Shaped Binaural For HS-103A	Two-in-one Tube System with Anatomically Shaped Binaural For HS-106A/107A/108A	Standard Binaural	Binaural for Rappaport

Accessory

Model No.	1711	1712A	1712B
			
	Colored PVC Tube	PVC Tube for Rappaport Type	PVC Transparent Tube

Accessory

Model No.	1741	1742	1743	1744
				
	Hard Plastic Ear Tip	Soft Plastic Ear Tip (Transparent)	Soft Plastic Ear Tip	Latex Ear Tip

Accessories

Accessory				
Model No.	1749	1746	1747	1748
				
	Translucent silicone Ear tip	Silicone Ear Tip	Deluxe Silicone Ear Tip	Antimicrobial Silicone Ear Tip
Accessory				
Model No.	1751	1752	1753	1754
				
	Plastic Membrane	Net Membrane	Transparent Membrane	Textured Membrane
Accessory				
Model No.	1755		1756	
				
	Mylar Membrane		Floating Diaphragm	
Accessory				
Model No.	174	173A	173B	
				
	Rectangular Name Tag	Metal Name Tag	Plastic Name Tag (for HS-101A)	
Accessory				
Model No.	176	177	178	
				
	Accessory Bag	Cannula for Rappaport tube	Accessories Box	



www.asian-medical.com

HONSUN (NANTONG) CO., LTD.

No.8 Tongxing Road, Economic & Technological Development
Area, 226009, Nantong City, China
TEL:0513-80580000 FAX:0513-80580080
E-mail: market@lordmed.com
www.asian-medical.com

SHANGHAI HONSUN INTERNATIONAL TRADE CO., LTD.

Rm. C. Floor 22. ZhiYuan Bldg., No. 768 Xietu
Road., Shanghai, P.R.China.
Tel: 86-21-63056696 Fax: 86-21-63036162
E-mail: market@lordmed.com
www.scian.com.cn